

Información general

Nombre de producto: ADAPTADOR DESTILACIÓN CON ENTRADA TERMÓMETRO

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Distancia "D"	Esmerilados Macho	Esmerilado Hembra
0100221411	200	14/23	14/23
0100221422	200	19/26	14/23
0100221433	260	24/29	14/23
0100221444	260	29/32	14/23

General Information

Product name: DISTILLING LINKS WITH THERMOMETER NECK

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Length "D"; (MM)	Cone Joints	Socket Joint
0100221411	200	14/23	14/23
0100221422	200	19/26	14/23
0100221433	260	24/29	14/23
0100221444	260	29/32	14/23