

Información general

Nombre de producto: ADAPTADOR TIPO CLAISEN DE 4 UNIONES

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Esmerilado P	Esmerilado M	Esmerilado N	Esmerilado O
0100222622	14/23	14/23	14/23	14/23
0100222641	29/32	19/26	19/26	14/23
0100222643	24/29	19/26	19/26	19/26
0100222644	29/32	29/32	29/32	14/23

Ficha de producto

Ref.	Esmerilado P	Esmerilado M	Esmerilado N	Esmerilado O
0100222645	24/29	24/29	24/29	24/29
0100222653	29/32	19/26	19/26	19/26
0100222655	29/32	29/32	29/32	29/32
0100222674	45/40	29/32	29/32	29/32



General Information

Product name: ADAPTER ACC. TO CLAISEN, 4 WAYS

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Joint P	Joint M	Joint N	Joint O
0100222622	14/23	14/23	14/23	14/23
0100222641	29/32	19/26	19/26	14/23
0100222643	24/29	19/26	19/26	19/26
0100222644	29/32	29/32	29/32	14/23

Product datasheet

Ref.	Joint P	Joint M	Joint N	Joint O
0100222645	24/29	24/29	24/29	24/29
0100222653	29/32	19/26	19/26	19/26
0100222655	29/32	29/32	29/32	29/32
0100222674	45/40	29/32	29/32	29/32

