

## Información general

Nombre de producto: COLUMNA LISA PARA RELLENO CON CAMISA DESMONTABLE

Fecha Rev.: Jul 11, 2023

Detalles:

## Fotografía



## Datos técnicos

Ref.	Longitud útil (MM)	Esmerilado Hembra/Macho
0100236021	200	14/23
0100236031	300	14/23
0100236034	300	29/32
0100236054	500	29/32
0100236084	800	29/32

## General Information

**Product name:** HEMPEL DISTILLATION COLUMN WITH REMOVABLE JACKET

**Rev. Date:** Jul 11, 2023

**Details:**

## Picture



## Technical Information

Ref.	Useful length (MM)	Socket/Cone Joints
0100236021	200	14/23
0100236031	300	14/23
0100236034	300	29/32
0100236054	500	29/32
0100236084	800	29/32