

Información general

Nombre de producto: MATRACES FONDO PLANO BOCA ESMERILADA

Fecha Rev.: Jul 11, 2023

Detalles:

Con boca esmerilada; en vidrio borosilicato.

Fotografía



Datos técnicos

| Ref. | Capacidad (ML) | Esmerilado |
|------------|----------------|------------|
| 0100111031 | 50 | 14/23 |
| 0100111032 | 50 | 19/26 |
| 0100111033 | 50 | 24/29 |

Ficha de producto

| Ref. | Capacidad (ML) | Esmerilado |
|------------|----------------|------------|
| 0100111034 | 50 | 29/32 |
| 0100111041 | 100 | 14/23 |
| 0100111042 | 100 | 19/26 |
| 0100111043 | 100 | 24/29 |
| 0100111044 | 100 | 29/32 |
| 0100111051 | 250 | 14/23 |
| 0100111052 | 250 | 19/26 |
| 0100111053 | 250 | 24/29 |
| 0100111054 | 250 | 29/32 |
| 0100111063 | 500 | 24/29 |
| 0100111064 | 500 | 29/32 |
| 0100111083 | 1000 | 24/29 |
| 0100111084 | 1000 | 29/32 |
| 0100111093 | 2000 | 24/29 |
| 0100111094 | 2000 | 29/32 |



General Information

Product name: FLAT BOTTOM FLASK WITH NS

Rev. Date: Jul 11, 2023

Details:

Short neck. Borosilicate glass.

Picture



Technical Information

| Ref. | Capacity (ML) | Joint |
|------------|---------------|-------|
| 0100111031 | 50 ML | 14/23 |
| 0100111032 | 50 ML | 19/26 |
| 0100111033 | 50 ML | 24/29 |

Product datasheet

| Ref. | Capacity (ML) | Joint |
|------------|---------------|-------|
| 0100111034 | 50 ML | 29/32 |
| 0100111041 | 100 ML | 14/23 |
| 0100111042 | 100 ML | 19/26 |
| 0100111043 | 100 ML | 24/29 |
| 0100111044 | 100 ML | 29/32 |
| 0100111051 | 250 ML | 14/23 |
| 0100111052 | 250 ML | 19/26 |
| 0100111053 | 250 ML | 24/29 |
| 0100111054 | 250 ML | 29/32 |
| 0100111063 | 500 | 24/29 |
| 0100111064 | 500 | 29/32 |
| 0100111083 | 1000 ML | 24/29 |
| 0100111084 | 1000 ML | 29/32 |
| 0100111093 | 2000 ML | 24/29 |
| 0100111094 | 2000 ML | 29/32 |

