


Información general

Nombre de producto: REFRIGERANTE DE DOBLE SUPERFICIE DIMROTH

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Longitud útil (MM)	Superf. intercambio (CM2)	Esmerilado Hembra/Macho
0100231534	300	850	29/32
0100231544	400	1135	29/32


General Information

Product name: COIL CONDENSER DOUBLE WALLED WITH OLIVES

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Useful length (MM)	Surface area (CM2)	Socket/Cone Joints
0100231534	300	850	29/32
0100231544	400	1135	29/32