

Información general

Nombre de producto: EMBUDO DE CARGA CON UNIÓN ESMERILADA INFERIOR, SIN GRADUAR

Fecha Rev.: Jul 11, 2023

Detalles:

Con unión esmerilada inferior, llave de vidrio o teflón y tapón de polietileno.

Fotografía



Datos técnicos

Ref.	Tapón	Tipo de llave	Capacidad (ML)	Esmerilado
0100238213	14/23	Vidrio	50	14/23
0100238214	19/26	Vidrio	100	14/23
0100238224	19/26	Vidrio	100	19/26

Ficha de producto

Ref.	Tapón	Tipo de llave	Capacidad (ML)	Esmerilado
0100238244	29/32	Vidrio	100	29/32
0100238245	29/32	Vidrio	250	29/32
0100238246	29/32	Vidrio	500	29/32
0100238248	29/32	Vidrio	1000	29/32
0100238414	19/26	PTFE cónica	100	14/23
0100238424	19/26	PTFE cónica	100	19/26
0100238444	29/32	PTFE cónica	100	29/32
0100238445	29/32	PTFE cónica	250	29/32
0100238446	29/32	PTFE cónica	500	29/32
0100238448	29/32	PTFE cónica	1000	29/32
0100238614	19/26	PTFE aguja	100	14/23
0100238624	19/26	PTFE aguja	100	19/26
0100238644	29/32	PTFE aguja	100	29/32
0100238645	29/32	PTFE aguja	250	29/32
0100238646	29/32	PTFE aguja	500	29/32
0100238648	29/32	PTFE aguja	1000	29/32



General Information

Product name: DROPPING FUNNEL, WITH BOTTOM JOINT, WITHOUT SCALE

Rev. Date: Jul 11, 2023

Details:

With bottom joint, stopcock and polyethylene stopper.

Picture



Technical Information

Ref.	Stopper	Stopcock	Capacity (ML)	Joint
0100238213	14/23	Glass	50 ML	14/23
0100238214	19/26	Glass	100 ML	14/23
0100238224	19/26	Glass	100 ML	19/26

Product datasheet

Ref.	Stopper	Stopcock	Capacity (ML)	Joint
0100238244	29/32	Glass	100 ML	29/32
0100238245	29/32	Glass	250 ML	29/32
0100238246	29/32	Glass	500	29/32
0100238248	29/32	Glass	1000 ML	29/32
0100238414	19/26	PTFE conical	100 ML	14/23
0100238424	19/26	PTFE conical	100 ML	19/26
0100238444	29/32	PTFE conical	100 ML	29/32
0100238445	29/32	PTFE conical	250 ML	29/32
0100238446	29/32	PTFE conical	500	29/32
0100238448	29/32	PTFE conical	1000 ML	29/32
0100238614	19/26	PTFE needle	100 ML	14/23
0100238624	19/26	PTFE needle	100 ML	19/26
0100238644	29/32	PTFE needle	100 ML	29/32
0100238645	29/32	PTFE needle	250 ML	29/32
0100238646	29/32	PTFE needle	500	29/32
0100238648	29/32	PTFE needle	1000 ML	29/32

