

Ficha de producto

Información general

Nombre de producto: EMBUDO DE CARGA CON UNIÓN ESMERILADA INFERIOR, GRADUADO

Fecha Rev.: Jul 11, 2023

Detalles:

Con unión esmerilada inferior, llave de vidrio o teflón y tapón de polietileno.

Fotografía



Datos técnicos

Ref.	Tapón	Tipo de llave	Capacidad (ML)	Esmerilado	Subdivisiones (ML)
0100238312	14/23	Vidrio	25	14/23	0.5
0100238314	19/26	Vidrio	100	14/23	2
0100238324	19/26	Vidrio	100	19/26	2

Ficha de producto

Ref.	Tapón	Tipo de llave	Capacidad (ML)	Esmerilado	Subdivisiones (ML)
0100238344	29/32	Vidrio	100	29/32	2
0100238345	29/32	Vidrio	250	29/32	5
0100238346	29/32	Vidrio	500	29/32	10
0100238348	29/32	Vidrio	1000	29/32	20
0100238514	19/26	PTFE cónica	100	14/23	2
0100238524	19/26	PTFE cónica	100	19/26	2
0100238544	29/32	PTFE cónica	100	29/32	2
0100238545	29/32	PTFE cónica	250	29/32	5
0100238546	29/32	PTFE cónica	500	29/32	10
0100238548	29/32	PTFE cónica	1000	29/32	20
0100238714	19/26	PTFE aguja	100	14/23	2
0100238724	19/26	PTFE aguja	100	19/26	2
0100238744	29/32	PTFE aguja	100	29/32	2
0100238745	29/32	PTFE aguja	250	29/32	5
0100238746	29/32	PTFE aguja	500	29/32	10
0100238748	29/32	PTFE aguja	1000	29/32	20

General Information

Product name: DROPPING FUNNEL, WITH BOTTOM JOINT, WITH SCALE

Rev. Date: Jul 11, 2023

Details:

With bottom joint, glass or PTFE stopcock and polyethylene stopper.

Picture



Technical Information

Ref.	Stopper	Stopcock	Capacity (ML)	Joint	Subdivisiones (ML)
0100238312	14/23	Glass	25 ML	14/23	0.5
0100238314	19/26	Glass	100 ML	14/23	2
0100238324	19/26	Glass	100 ML	19/26	2

Product datasheet

Ref.	Stopper	Stopcock	Capacity (ML)	Joint	Subdivisiones (ML)
0100238344	29/32	Glass	100 ML	29/32	2
0100238345	29/32	Glass	250 ML	29/32	5
0100238346	29/32	Glass	500	29/32	10
0100238348	29/32	Glass	1000 ML	29/32	20
0100238514	19/26	PTFE conical	100 ML	14/23	2
0100238524	19/26	PTFE conical	100 ML	19/26	2
0100238544	29/32	PTFE conical	100 ML	29/32	2
0100238545	29/32	PTFE conical	250 ML	29/32	5
0100238546	29/32	PTFE conical	500	29/32	10
0100238548	29/32	PTFE conical	1000 ML	29/32	20
0100238714	19/26	PTFE needle	100 ML	14/23	2
0100238724	19/26	PTFE needle	100 ML	19/26	2
0100238744	29/32	PTFE needle	100 ML	29/32	2
0100238745	29/32	PTFE needle	250 ML	29/32	5
0100238746	29/32	PTFE needle	500	29/32	10
0100238748	29/32	PTFE needle	1000 ML	29/32	20

