

Información general

Nombre de producto: EMBUDO DE DECANTACIÓN CÓNICO, SIN GRADUAR

Fecha Rev.: Jul 11, 2023

Detalles:

Con llave y tapón de polietileno.

Fotografía



Datos técnicos

Ref.	Tapón	Tipo de llave	Capacidad (ML)	Ø Exterior vástago (MM)
0100239403	14/23	Vidrio	50	8
0100239404	19/26	Vidrio	100	9
0100239405	29/32	Vidrio	250	9

Ficha de producto

Ref.	Tapón	Tipo de llave	Capacidad (ML)	Ø Exterior vástago (MM)
0100239406	29/32	Vidrio	500	10
0100239408	29/32	Vidrio	1000	10
0100239409	29/32	Vidrio	2000	12
0100239502	14/23	PTFE cónica	25	8
0100239503	14/23	PTFE cónica	50	8
0100239504	19/26	PTFE cónica	100	9
0100239505	29/32	PTFE cónica	250	9
0100239506	29/32	PTFE cónica	500	10
0100239508	29/32	PTFE cónica	1000	10
0100239509	29/32	PTFE cónica	2000	12
0100239511	45/40	PTFE cónica	5000	10
0100239603	14/23	PTFE aguja	50	8
0100239604	19/26	PTFE aguja	100	9
0100239605	29/32	PTFE aguja	250	9
0100239606	29/32	PTFE aguja	500	10
0100239608	29/32	PTFE aguja	1000	10
0100239609	29/32	PTFE aguja	2000	12



General Information

Product name: SEPARATING FUNNEL, WITHOUT SCALE

Rev. Date: Jul 11, 2023

Details:

With stopcock and polyethylene stopper.

Picture



Technical Information

Ref.	Stopper	Stopcock	Capacity (ML)	OD Stem (MM)
0100239403	14/23	Glass	50 ML	8
0100239404	19/26	Glass	100 ML	9
0100239405	29/32	Glass	250 ML	9

Product datasheet

Ref.	Stopper	Stopcock	Capacity (ML)	OD Stem (MM)
0100239406	29/32	Glass	500	10
0100239408	29/32	Glass	1000 ML	10
0100239409	29/32	Glass	2000 ML	12
0100239502	14/23	PTFE conical	25 ML	8
0100239503	14/23	PTFE conical	50 ML	8
0100239504	19/26	PTFE conical	100 ML	9
0100239505	29/32	PTFE conical	250 ML	9
0100239506	29/32	PTFE conical	500	10
0100239508	29/32	PTFE conical	1000 ML	10
0100239509	29/32	PTFE conical	2000 ML	12
0100239511	45/40	PTFE conical	5000	10
0100239603	14/23	PTFE needle	50 ML	8
0100239604	19/26	PTFE needle	100 ML	9
0100239605	29/32	PTFE needle	250 ML	9
0100239606	29/32	PTFE needle	500	10
0100239608	29/32	PTFE needle	1000 ML	10
0100239609	29/32	PTFE needle	2000 ML	12

