

Información general

Nombre de producto: FRASCOS DE BOCA ANCHA

Fecha Rev.: Jul 11, 2023

Detalles:

Frasco cilíndrico en polietileno de color natural. Con los orificios del tapón y la anilla del cuerpo, permite sellar el frasco mediante alambre. Los frascos se suministran sin roscar. Cuerpo de PE y tapa de PP.

Fotografía



Datos técnicos

Ref.	Ø Boca (MM)	Capacidad (ML)	Unidades (Pack)
10FRTE0060	24	50	10
10FRTE0125	24	100	10
10FRTE0250	34.5	250	10
10FRTE0500	34.5	500	10

Ficha de producto

Ref.	Ø Boca (MM)	Capacidad (ML)	Unidades (Pack)
10FRTE1000	45	1000	10
10FRTE2000	45	2000	10



General Information

Product name: FLASKS WITH SCREW CAP, WIDE MOUTH

Rev. Date: Jul 11, 2023

Details:

Cylindrical bottle made of high density polyethylene, natural color. You can seal the bottle by a wire with the holes of the stopper and the o-ring of the body. The bottles are supplied without closing. Body (PE), cap (PP)

Picture



Technical Information

Ref.	Ø Mouth (MM)	Capacity (ML)	Units (Pack)
10FRTE0060	24	50 ML	10
10FRTE0125	24	100 ML	10
10FRTE0250	34.5	250 ML	10
10FRTE0500	34.5	500	10

Product datasheet

Ref.	Ø Mouth (MM)	Capacity (ML)	Units (Pack)
10FRTE1000	45	1000 ML	10
10FRTE2000	45	2000 ML	10

