

## Información general

Nombre de producto: REACTORES CILÍNDRICOS FONDO REDONDO CON CAMISA Y LLAVE

Fecha Rev.: Jul 11, 2023

Detalles:

## Fotografía



## Datos técnicos

Ref.	Ø Ext. Reactor (mm)	Ø Int. Reactor (mm)	Llave (mm)	Conector	Capacidad (ML)	Altura Total (mm)	DN Boca
0100802400	90	60	10	GL14	100	245	60
0100802402	90	60	10	GL14	250	300	60
0100802404	90	60	10	GL14	500	385	60

## Ficha de producto

Ref.	Ø Ext. Reactor (mm)	Ø Int. Reactor (mm)	Llave (mm)	Conector	Capacidad (ML)	Altura Total (mm)	DN Boca
0100802406	90	60	10	GL14	1000	385	60
0100802416	130	99	10	GL18	250	275	100
0100802418	130	99	10	GL18	500	295	100
0100802420	130	99	10	GL18	1000	360	100
0100802422	130	99	10	GL18	2000	490	100
0100802424	190	150	10	GL18	3000	440	100
0100802426	190	150	10	GL18	4000	495	100
0100802440	190	150	10	GL18	2000	375	150
0100802442	190	150	10	GL18	3000	430	150
0100802444	190	150	10	GL18	4000	485	150
0100802446	190	150	10	GL18	5000	540	150
0100802448	190	150	10	GL18	6000	595	150
0100802450	315	240	20	GL18	10000	640	150
0100802466	250	200	10	KF15	6000	535	200
0100802468	250	200	20	KF15	10000	665	200

# Product datasheet

## General Information

**Product name:** CYLINDRICAL REACTION VESSELS, ROUND BOTTOM, WITH JACKET AND STOPCOCK

**Rev. Date:** Jul 11, 2023

**Details:**

## Picture



## Technical Information

Ref.	Vessel Outer Ø (mm)	Vessel Inner Ø (mm)	Valve (mm)	Connector	Capacity (ML)	Total Height (mm)	Laboratory Flange DN
0100802400	90	60	10	GL14	100 ML	245	60
0100802402	90	60	10	GL14	250 ML	300	60
0100802404	90	60	10	GL14	500	385	60

## Product datasheet

Ref.	Vessel Outer $\varnothing$ (mm)	Vessel Inner $\varnothing$ (mm)	Valve (mm)	Connector	Capacity (ML)	Total Height (mm)	Laboratory Flange DN
0100802406	90	60	10	GL14	1000 ML	385	60
0100802416	130	99	10	GL18	250 ML	275	100
0100802418	130	99	10	GL18	500	295	100
0100802420	130	99	10	GL18	1000 ML	360	100
0100802422	130	99	10	GL18	2000 ML	490	100
0100802424	190	150	10	GL18	3000 ML	440	100
0100802426	190	150	10	GL18	4000	495	100
0100802440	190	150	10	GL18	2000 ML	375	150
0100802442	190	150	10	GL18	3000 ML	430	150
0100802444	190	150	10	GL18	4000	485	150
0100802446	190	150	10	GL18	5000	540	150
0100802448	190	150	10	GL18	6000	595	150
0100802450	315	240	20	GL18	10000	640	150
0100802466	250	200	10	KF15	6000	535	200
0100802468	250	200	20	KF15	10000	665	200