

## Información general


**Nombre de producto:** EJE AGITADOR DE PALAS DE HÉLICE

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

- Químicamente inerte.
- Fácil de limpiar.
- Evita la contaminación cruzada.
- Resistente a roturas y arañazos, y hasta 280 °C.
- Los modelos de eje expuesto (EX) tienen un núcleo de acero inoxidable, con un extremo expuesto de longitud de 50 mm para una fijación más rígida en el mandril.

## Fotografía



## Datos técnicos

## Ficha de producto

Ref.	Ø Eje (mm)	Ø Rotor (mm)	Longitud (mm)
10AGHE0630	6	50	500
10AGHE0835	8	50	500
10AGHE1010	10	70	550
10AGHE1215	12	80	650EX
10AGHE1220	16	100	750EX



## General Information


**Product name:** SHAFT STIRRER, SCREW PROPELLER

**Rev. Date:** Jul 11, 2023

**Details:**

- Chemically inert.
- Easy to clean.
- No cross contamination.
- Resistant to breakage, scratches and heat up to 280 °C.
- Exposed stainless steel core models (EX) provide extra rigid clamping in the drive chuck, with an exposed length of 50 mm.

## Picture



## Technical Information

## Product datasheet

Ref.	Shaft $\varnothing$ (mm)	Rotor $\varnothing$ (mm)	Length (mm)
10AGHE0630	6	50	500
10AGHE0835	8	50	500
10AGHE1010	10	70	550
10AGHE1215	12	80	650EX
10AGHE1220	16	100	750EX

