

Información general

Nombre de producto: PICNÓMETRO GAY-LUSSAC CON TAPÓN

Fecha Rev.: Jul 11, 2023

Detalles:

Con tapón capilar esmerilado 10/19, con o sin calibración.

Fotografía



Datos técnicos

Ref.	Capacidad (ML)	Calibración
0100521400	5	Sin calibrar
0100521401	10	Sin calibrar
0100521402	25	Sin calibrar

Ficha de producto

Ref.	Capacidad (ML)	Calibración
0100521403	50	Sin calibrar
0100521404	100	Sin calibrar
0100521500	5	Calibrado con certificado individual
0100521501	10	Calibrado con certificado individual
0100521502	25	Calibrado con certificado individual
0100521503	50	Calibrado con certificado individual
0100521504	100	Calibrado con certificado individual



General Information

Product name: PYCNOMETERS GAY-LUSSAC, WITH STOPPER

Rev. Date: Jul 11, 2023

Details:

With capillary stopper NS 10/19, with or without calibration.

Picture



Technical Information

Ref.	Capacity (ML)	Calibration
0100521400	5 ML	Unadjusted
0100521401	10 ML	Unadjusted
0100521402	25 ML	Unadjusted

Product datasheet

Ref.	Capacity (ML)	Calibration
0100521403	50 ML	Unadjusted
0100521404	100 ML	Unadjusted
0100521500	5 ML	Adjusted with single certificate
0100521501	10 ML	Adjusted with single certificate
0100521502	25 ML	Adjusted with single certificate
0100521503	50 ML	Adjusted with single certificate
0100521504	100 ML	Adjusted with single certificate

